


The RMS Quick Release Multi-Mounting Frame



ASSESSMENT / ORDER FORM

(QMF 49) Mar. 2006 

MAIN FEATURES:

The Quick Release Multi-Mounting Frame system was developed by RMS to offer singular or combined use for:-

- Headrest attachment - Shoulder Strap mounting – Backrest Tensioning.
- Fully modular, with height and width adjustment, allowing easy positioning to suit individual requirements.
- Quick Release mounting system enables easy removal for folding of wheelchair.
- Frame finished in durable chromium plating.
- Ideally suited for users with two wheelchairs, as with this system, it is only necessary to purchase and fit additional mounting brackets (MF2#, see parts list below) and the main cross frame with headrest for example, can be transferred quit easily between wheelchairs.
- In October 1991 the QRMMF, complete with a conventional Headrest, was attached to a generic style of wheelchair that was successfully “Crash Tested” to the frontal impact requirements of ISO 7176/19 at Millbrook Proving Ground.



PARTS LIST:

Mounting Frame (width adjustable)		Tick
• Suitable for wheelchair width 10” to 14” (255 to 355mm)	MF1/10
• Suitable for wheelchair width 15” to 20” (380 to 500mm)	MF1/15
Backrest Mounting Clamps – pair (to fit 7/8” – 22mm backrest tube)		
• 80mm above Backrest Canvas top edge	MF2/80
• 120mm above Backrest Canvas top edge	MF2/120
• 160mm above Backrest Canvas top edge	MF2/160
• For attachment to Backrests with wrap-around canvas	MF2/C
Reduction Spacers - pair		
• To reduce 7/8” to 3/4” (22mm –19mm)	MF3
Shoulder Harness Mountings – pair		
• Harness positioning collars	MF4
Note: Spacer Shim MF5 supplied with MF4 Kit.		
Webbing Guides – pair		
• 1” (25mm)	MP001
• 2” (40-50mm)	MP003
Headrest Mounting		
• Receiver socket for headrest vertical stem	MF6

WHEELCHAIR SERVICE:	CLIENT NAME / REF:	Need some Help? Call our
CONTACT NAME:	DATE:	Technical Help-line on 01634-572513