

KAYE

POSTURE CONTROL

WALKER



KAYE WALKER FEATURES

- Designed to make walking less energy-consuming
- Can be included in rehabilitation programs, for both indoor or outdoor activities
- Fully height adjustable frame with multiple positioning options
- Fixed, activity and directional locking swivel wheels options available
- Available in seven different sizes including a heavy duty frame option
- A range of accessories available to best meet the unique needs of every user
- Easy folding, lightweight and easy to carry

BENEFITS OF THE KAYE WALKER

The Kaye Posture Control Walker is the perfect product to be used as part of a walking development programme in schools, home, or out in the community.

A product used by the NHS and equipment providers for over 20 years, the Kaye Walker offers a wide range of parts and accessories to provide users with a stable walking platform to work from.

The stainless steel and aluminium frame is sturdy and secure, yet remains extremely light weight; allowing for easy folding to go into the boot of the car; but strong enough to withstand years of constant use.



WHY THIS DESIGN?

The Kaye Posture Control Walker is designed to make walking less energy consuming and maximise the potential for walking. The design and function of the walker encourage hip and trunk extension to develop muscle tone and balance necessary for walking, whilst offering independence and more sociability to users.



TECHNICAL SPECIFICATIONS

Helps to improve the users natural bodily functions.

Promotes users independence, social interaction and engagement.

Encourages the maintenance of a vertical postural alignment.

Contributes to strengthening upper and lower body muscles.

	WALKER HEIGHT STANDARD WHEEL	WALKER HEIGHT ACTIVITY WHEEL	HAND GRIP WIDTH	FRAME WIDTH	FRAME LENGTH	HEIGHT OF USER	MAX USER WEIGHT	FRAME WEIGHT
W1/2	37 - 46	N/A	34	58 - 60	52 - 59	up to 95	27 kg	3.00 kg
W1	41 - 55	45 - 59	34	60 - 62	56 - 62	91 - 122	27 kg	3.00 kg
W2	48 - 64	52 - 68	34	58 - 60	69 - 83	107 - 137	39 kg	3.85 kg
W3	59 - 78	63 - 82	38	62 - 64	77 - 89	129 - 152	60 kg	4.85 kg
W4L	72 - 92	76 - 96	40	65 - 67	84 - 97	150 - 179	82 kg	5.50 kg
W4H	81 - 100	N/A	40	67 - 69	90 - 102	150 - 185	95 kg	7.00 kg
W5	90 - 105	94 - 109	47	65 - 67	102 - 108	175 - 192	113 kg	9.75 kg



Thompson House, Unit 10, Styles Close,
Sittingbourne, Kent, ME10 3BF.

Tel: +44 (0) 1795 477280

Fax: +44 (0) 1795 229692

Email: sales@rms-rehab.co.uk

www.rms-rehab.co.uk

